



PRODUCT SPECIFICATIONS

Product Name:	MAb to HRV	
Description:	Monoclonal Antibody to Human Rotavirus	
Catalog#:	M012103M	Lot #: MY21M1103
Isotype:	IgG2a	
Immunogen:	Purified Rotavirus	
Source:	Mouse ascites	
Format:	Purified, Liquid	
Concentration :	3.42 mg/mL (+/- 10 %) (OD280nm, E ^{0.1%} = 1.4)	
Purity :	≥90% (SDS-PAGE)	
Product Buffer:	PBS, PH7.4	
Preservative:	0.05% Proclin 300	
Applications:	Suitable for use in immunoassays. Each laboratory should determine an optimum working titer for use in its particular applications.	
Storage:	Store at 2-8°C.	
Production Date:	03/2021	
Expiry Date:	03/2024	

Note: Centrifuge before opening to ensure complete recovery of vial contents

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Web: www.firsunbioscience.com

Add: Building C5, Hongfeng Science and Technology Park, Qixia District, Nanjing, Jiangsu Province, China, 210000

Tel:86-025-83108003 Email: order@firsunbioscience.com