



PRODUCT SPECIFICATIONS

Product Name:	Influenza B-N Protein	
Description:	Recombinant Influenza B nucleocapsid protein with N-terminus HIS tag.	
Catalog #:	A011301R	Lot #: AY22M1905
Source:	<i>Escherichia coli</i>	
Format:	Purified, Liquid	
Concentration:	2.75mg/mL (+/- 10 %)	
Purity:	≥85% (SDS-PAGE)	
Product Buffer:	PBS, PH7.4	
Preservative:	0.05%Proclin 300	
Applications:	Suitable for use in immunoassays. Each laboratory should determine an optimum working titer for use in its particular applications.	
Inactivation:	Whole infective virus is inactivated with formalin. Confirm by attempting growth under original culture conditions.	
Storage:	Store at 2-8°C for two weeks and -20°C for long time. Avoid multiple freeze/thaw cycles.	
production date:	05/2022	
Expiry Date:	05/2025	

Note: Centrifuge before opening to ensure complete recovery of vial contents

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Web: www.firsunbioscience.com

Add: Building C5, Hongfeng Science and Technology Park, Qixia District, Nanjing, Jiangsu Province, China, 210000

Tel:86-025-83108003 Email: order@firsunbioscience.com