



PRODUCT SPECIFICATIONS

Product Name:	MAb to nCoV-N	
Description:	Monoclonal Antibody to novel coronavirus nucleocapsid protein	
Catalog #:	M011503M	Lot #: MY22M1702
Isotype:	IgG1	
Immunogen:	Purified nCoV-N protein	
Source:	cell culture	
Format:	Purified, Liquid	
Concentration:	3.56 mg/mL (+/- 10 %) (OD280nm, E ^{0.1%} = 1.4)	
Purity:	≥90% (SDS-PAGE)	
Product Buffer:	PBS, PH7.4	
Preservative:	0.05% ProClin300	
Applications:	Suitable for use in immunoassays. Each laboratory should determine an optimum working titer for use in its particular applications.	

Capture

M011503M

Detection

M011504M

Storage:	Store at 2-8°C for two weeks and -20°C for long time. Avoid multiple freeze/thaw cycles.
production date :	02/2022
Expiry Date:	02/2025

Note: Centrifuge before opening to ensure complete recovery of vial contents

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Web: www.firsunbioscience.com

Add: Building C5, Hongfeng Science and Technology Park, Qixia District, Nanjing, Jiangsu Province, China, 210000

Tel:86-025-83108003 Email: order@firsunbioscience.com