



PRODUCT SPECIFICATIONS

Product Name:	MAb to α -IFN	
Description:	Monoclonal Antibody to Interferon alpha	
Catalog #:	M031002M	Lot #: [Lot specific]
Isotype:	IgG1	
Immunogen:	Purified α -IFN protein	
Source:	Mouse ascites	
Format:	Purified, Liquid	
Concentration :	[Lot specific]mg/mL (+/- 10 %) (OD280nm, $E^{0.1\%} = 1.4$)	
Purity :	$\geq 90\%$ (SDS-PAGE)	
Product Buffer:	PBS, PH7.4	
Preservative:	0.05% Proclin 300	
Applications:	Suitable for use in immunoassays. Each laboratory should determine an optimum working titer for use in its particular applications.	

	Capture	Detection
	M031002M	M031001M
Storage:	Store at 2-8°C for two weeks and -20°C for long time. Avoid multiple freeze/thaw cycles.	
Production Date:	[Lot specific]	
Expiry Date:	[Lot specific]	

Note: Centrifuge before opening to ensure complete recovery of vial contents

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Web: www.firsunbioscience.com

Add: Building C5, Hongfeng Science and Technology Park, Qixia District, Nanjing, Jiangsu Province, China, 210000

Tel:86-025-83108003 Email: order@firsunbioscience.com