



PRODUCT SPECIFICATIONS

Product Name:	MAb to CEA		
Description:	Monoclonal Antibody to Carcinoembryonic antigen		
Catalog #:	M021802M	Lot #:	MY21M1210
Isotype:	IgG1		
Immunogen:	Purified CEA protein		
Source:	Mouse ascites		
Format:	Purified, Liquid		
Concentration :	3.2 mg/mL (+/- 10 %) (OD280nm, E ^{0.1%} = 1.4)		
Purity :	≥90% (SDS-PAGE)		
Product Buffer:	PBS, PH7.4		
Preservative:	0.05% Proclin 300		
Applications:	Suitable for use in immunoassays. Each laboratory should determine an optimum working titer for use in its particular applications.		

Capture

Detection

M021802M

M021803M/M021804M/M021805M

/M021806M

Storage:	Store at 2-8°C for two weeks and -20°C for long time. Avoid multiple freeze/thaw cycles.
Production Date:	10/2021
Expiry Date:	10/2024

Note: Centrifuge before opening to ensure complete recovery of vial contents

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Web: www.firsunbioscience.com

Add: Building C5, Hongfeng Science and Technology Park, Qixia District, Nanjing, Jiangsu Province, China, 210000

Tel:86-025-83108003 Email: order@firsunbioscience.com