



PRODUCT SPECIFICATIONS

Product Name:	HPV-16 Protein	
Description:	Recombinant Human Papilloma Virus Type 16 (HPV16), L1 Protein with C-terminal HIS tag	
Catalog #:	A020802R	Lot #: AY21M3112
Source:	<i>Escherichia coli</i>	
Format:	Purified, Liquid	
Concentration :	1.0 mg/mL (+/- 10 %)	
Purity :	≥90% (SDS-PAGE)	
Product Buffer:	50mM Tris-HCl , 100mM NaCl , 5% glycine , 10% glycerol, pH6.5	
Preservative:	None	
Applications:	Suitable for use in immunoassays. Each laboratory should determine an optimum working titer for use in its particular applications.	
Storage:	Store at 2-8°C for two weeks and -20°C for long time. Avoid multiple freeze/thaw cycles.	
Production Date:	12/2021	
Expiry Date:	12/2024	

Note: Centrifuge before opening to ensure complete recovery of vial contents

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Web: www.firsunbioscience.com

Add: Building C5, Hongfeng Science and Technology Park, Qixia District, Nanjing, Jiangsu Province, China, 210000

Tel:86-025-83108003 Email: order@firsunbioscience.com