




PRODUCT SPECIFICATIONS

Product Name:	HPV16 E7 Protein	
Description:	Recombinant Human Papilloma Virus Type 16 (HPV16), E7 Protein with C-terminal HIS tag	
Catalog #:	A020804R	Lot #: AY22M0211
Source:	<i>Escherichia coli</i>	
Format:	Purified, Liquid	
Concentration:	1.0 mg/mL (+/- 10 %)	
Purity:	≥90% (SDS-PAGE)	
Product Buffer:	Tris ,PH8.5	
Preservative:	None	
Applications:	Suitable for use in immunoassays. Each laboratory should determine an optimum working titer for use in its particular applications.	
Storage:	Store at 2-8°C for two weeks and -20°C for long time. Avoid multiple freeze/thaw cycles.	
Production Date:	11/2022	
Expiry Date:	11/2025	
SDS-PAGE:		

Note: Centrifuge before opening to ensure complete recovery of vial contents

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Web: www.firsunbioscience.com

Add: Building C5, Hongfeng Science and Technology Park, Qixia District, Nanjing, Jiangsu Province, China, 210000

Tel:86-025-83108003 Email: order@firsunbioscience.com