



## PRODUCT SPECIFICATIONS

<b>Product Name:</b>	MAb to C-IV	
<b>Description:</b>	Monoclonal Antibody to Collage Type IV	
<b>Catalog #:</b>	M080902M	<b>Lot #:</b> [Lot specific]
<b>Isotype:</b>	IgG1	
<b>Immunogen:</b>	Purified C-IV protein	
<b>Source:</b>	Mouse ascites	
<b>Format:</b>	Purified, Liquid	
<b>Concentration :</b>	[Lot specific] mg/mL (+/- 10 %) (OD280nm, E <sup>0.1%</sup> = 1.4)	
<b>Purity :</b>	≥90% (SDS-PAGE)	
<b>Product Buffer:</b>	PBS, PH7.4	
<b>Preservative:</b>	0.05% Proclin 300	
<b>Applications:</b>	Suitable for use in immunoassays. Each laboratory should determine an optimum working titer for use in its particular applications.	
	<b>Capture</b>	<b>Detection</b>
	M080902M	M080901M
<b>Storage:</b>	Store at 2-8°C for two weeks and -20°C for long time. Avoid multiple freeze/thaw cycles.	
<b>Production Date:</b>	[Lot specific]	
<b>Expiry Date:</b>	[Lot specific]	

**Note:** Centrifuge before opening to ensure complete recovery of vial contents

### FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Web: [www.firsunbioscience.com](http://www.firsunbioscience.com)

Add: Building C5, Hongfeng Science and Technology Park, Qixia District, Nanjing, Jiangsu Province, China, 210000

Tel:86-025-83108003      Email: [order@firsunbioscience.com](mailto:order@firsunbioscience.com)