



PRODUCT SPECIFICATIONS

Product Name:	CG-BSA Conjugate	
Description:	Cholyglycine-BSA Conjugate	
Catalog #:	A080102R	Lot #: AY21M1012
Source:	Synthesis	
Format:	Purified, Liquid	
Concentration :	20.02 mg /mL (+/- 10 %)	
Purity :	None	
Product Buffer:	PBS, PH7.4	
Preservative:	0.05% Proclin 300	
Applications:	Suitable for use in immunoassays. Each laboratory should determine an optimum working titer for use in its particular applications.	
Storage:	Store at 2-8°C in dark.	
Production Date:	12/2021	
Expiry Date:	12/2024	

Note: Centrifuge before opening to ensure complete recovery of vial contents

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Web: www.firsunbioscience.com

Add: Building C5, Hongfeng Science and Technology Park, Qixia District, Nanjing, Jiangsu Province, China, 210000

Tel:86-025-83108003 Email: order@firsunbioscience.com