



PRODUCT SPECIFICATIONS

Product Name:	MAb to FCV		
Description:	Monoclonal Antibody to Feline Calicivirus		
Catalog #:	M100703M	Lot #:	MY22M2908
Isotype:	IgG1		
Immunogen:	Purified FCV Protein		
Source:	Mouse ascites		
Format:	Purified, Liquid		
Concentration:	2.1 mg/mL (+/- 10 %) (OD280nm, $E^{0.1\%} = 1.4$)		
Purity:	≥90% (SDS-PAGE)		
Product Buffer:	PBS, PH7.4		
Preservative:	0.05% Proclin 300		
Applications:	Suitable for use in immunoassays. Each laboratory should determine an optimum working titer for use in its particular applications.		
Storage:	Store at 2-8°C for two weeks and -20°C for long time. Avoid multiple freeze/thaw cycles.		
production date:	08/2022		
Expiry Date:	08/2025		

Note: Centrifuge before opening to ensure complete recovery of vial contents

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Web: www.firsunbioscience.com

Add: Building C5, Hongfeng Science and Technology Park, Qixia District, Nanjing, Jiangsu Province, China, 210000

Tel:86-025-83108003 Email: order@firsunbioscience.com