



PRODUCT SPECIFICATIONS

Product Name:	MAb to FCV	
Description:	Monoclonal Antibody to Feline Calicivirus	
Catalog #:	M100707M	Lot #: MY22M2111
Isotype:	IgG1	
Immunogen:	Purified FCV protein	
Source:	Mouse ascites	
Format:	Purified, Liquid	
Concentration:	3.0 mg/mL (+/- 10 %) (OD280nm, E ^{0.1%} = 1.4)	
Purity:	≥90% (SDS-PAGE)	
Product Buffer:	PBS, PH7.4	
Preservative:	0.05% Proclin 300	
Applications:	Suitable for use in immunoassays. Each laboratory should determine an optimum working titer for use in its particular applications.	
	Capture	Detection
	M100705M/M100708M	M100707M
Storage:	Store at 2-8°C for two weeks and -20°C for long time. Avoid multiple freeze/thaw cycles.	
production date:	11/2022	
Expiry Date:	11/2025	

Note: Centrifuge before opening to ensure complete recovery of vial contents

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Web: www.firsunbioscience.com

Add: Building C5, Hongfeng Science and Technology Park, Qixia District, Nanjing, Jiangsu Province, China, 210000

Tel:86-025-83108003 Email: order@firsunbioscience.com