




PRODUCT SPECIFICATIONS

| | | |
|-------------------------|---|-------------------------|
| Product Name: | FPV Protein | |
| Description: | Feline Panleukopenia Virus Protein with C-terminal HIS tag | |
| Catalog #: | A100603R | Lot #: AY22M3108 |
| Source: | <i>Escherichia coli</i> | |
| Format: | Purified, Liquid | |
| Concentration: | 1.0 mg/mL (+/- 10 %) | |
| Purity: | ≥90% (SDS-PAGE) | |
| Product Buffer: | Tris , PH8.5 | |
| Preservative: | 0.05% Proclin 300 | |
| Applications: | Suitable for use in immunoassays. Each laboratory should determine an optimum working titer for use in its particular applications. | |
| Storage: | Store at 2-8°C for two weeks and -20°C for long time . Avoid multiple freeze/thaw cycles. | |
| production date: | 08/2022 | |
| Expiry Date: | 08/2025 | |

SDS-PAGE:



Web: www.firsunbioscience.com

Add: Building C5, Hongfeng Science and Technology Park, Qixia District, Nanjing, Jiangsu Province, China, 210000

Tel:86-025-83108003 Email: order@firsunbioscience.com



Note: Centrifuge before opening to ensure complete recovery of vial contents

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Web: www.firsunbioscience.com

Add: Building C5, Hongfeng Science and Technology Park, Qixia District, Nanjing, Jiangsu Province, China, 210000

Tel:86-025-83108003 Email: order@firsunbioscience.com