



PRODUCT SPECIFICATIONS

Product Name:	CAV-1 Protein	
Description:	Canine Adenovirus -1 Protein	
Catalog#:	A101202N	Lot #: AY22M1508
Source:	<i>F81</i>	
Format:	Purified, Liquid	
Concentration :	13 mg/mL (+/- 10 %)	
Purity :	None	
Product Buffer:	PBS, PH7.4	
Preservative:	N/A	
Applications:	Suitable for use in immunoassays. Each laboratory should determine an optimum working titer for use in its particular application.	
Storage:	Store at 2-8°C for two weeks and -20°C for long time. Avoid multiple freeze/thaw cycles.	
Production Date:	08/2022	
Expiry Date:	08/2025	

Note: Centrifuge before opening to ensure complete recovery of vial contents

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Web: www.firsunbioscience.com

Add: Building C5, Hongfeng Science and Technology Park, Qixia District, Nanjing, Jiangsu Province, China, 210000

Tel:86-025-83108003 Email: order@firsunbioscience.com