



PRODUCT SPECIFICATIONS

Product Name:	PAb to Rabbit anti CCRP	
Description:	Polyclonal Antibody Rabbit to Canine C-Reactive Protein	
Catalog #:	P100501R	Lot #: PY23M2505
Isotype:	IgG	
Immunogen:	Purified CCRP protein	
Source:	Rabbit Serum	
Format:	Purified, Liquid	
Concentration:	14.0 mg/mL (+/- 10 %) (OD280nm, E ^{0.1%} = 1.4)	
Purity:	≥90% (SDS-PAGE)	
Product Buffer:	PBS, pH 7.4	
Preservative:	0.05% Proclin 300	
Applications:	Suitable for use in immunoassays. Each laboratory should determine an optimum working titer for use in its particular applications. Recommended pairs for immunoassays:	
Storage:	Store at 2-8°C for two week and -20°C for long time. Avoid multiple freeze/thaw cycles.	
Production Date:	05/2023	
Expiry Date:	05/2026	

Note: Centrifuge before opening to ensure complete recovery of vial contents

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Web: www.firsunbioscience.com

Add: Building C5, Hongfeng Science and Technology Park, Qixia District, Nanjing, Jiangsu Province, China, 210000

Tel:86-025-83108003 Email: order@firsunbioscience.com